Click to edit Master

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SIMPLIFY YOUR STARTUP
FOR TARGETED VETERINARY DRUG SCREENING
QQQ LC/MS Veterinary Drug Application Kit
 The Measure of Confidence

Anilent LC/MS Application Kits

Minimize the need for tedious manual method development with pretested methods and a triggered MRM (tMRM) database and library

Analyzing veterinary drug compounds is challenging for two reasons: low concentrations in difficult matrices and a large number of analytes to be monitored and quantified. With so many variables to consider, it can be difficult to find a solid starting point for method development.

A faster, easier way to develop customized screening methods

Agilent's Veterinary Drug Screening tMRM Application Kit is truly unique, because much of the development work has already been completed. The kit features easy-to-use examples that show you how to set up screening methods and quickly adapt them to your specific needs. It also includes:

- A tMRM database and library for several hundred Veterinary Drug compounds that includes compound names, up to 10 MRM transitions, fragmentor voltages, collision energies, and retention times for each database compound for reliable forensic screening with tMRM library verification.

The following components are included – saving you time and money:

- Agilent's Veterinary Drug tMRM database and library
- Comprehensive standard test mix with more than 140 compounds
- Agilent LC columns ideally suited for Rapid Resolution Liquid Chromatography (RRLC) and Ultra High Pressure Liquid Chromatography (UHPLC) methods
- Quick-start guide and Application Note that show you how to run the test mixes and create tMRM methods
- CD-ROM with examples of screening methods, data files, and reports
- On site application-based training to assist with method setup for quick and effective startup


Ordering information:

Agilent's QQQ LC/MS Veterinary Drug Application Kit (G1735AA)

The following are required but not included with the G1735AA kit:

- Agilent 1200 Series RRLC System or Agilent 1290 Infinity LC
- Agilent 6400 Series Triple Quadrupole LC/MS system
- Agilent MassHunter Acquisition Software B.06 or higher

Pre-developed examples help you implement your screening method in a fraction of the time



tMRM acquisition produces a tMRM product ion spectrum with up to 10 spectral peaks per compound present. These spectra are searched against a tMRM reference library for compound verification. A typical tMRM spectrum is shown here.


Run segment																	
Compound #/mol	Compound Name	RT/yr	Pressure (torr)	MS1 Res	Product Ion	MS1 Res	Pressure	Trigger	Threshold	Ref Time (s)	Exp Time (s)	Fragment	Cation Charge	Collisional Energy	Collisional Voltage	Isobar	Signal
M	Dimethylsiloxane	1.0	405.1 Torr	441.1 Torr	SP	SP	SP	SP	500	1.4	1	125	12	7.7kV	0	0	0
	Dimethylsiloxane	1.0	405.1 Torr	441.1 Torr	SP	SP	SP	SP	500	1.4	1	125	20	7.7kV	0	0	0
	Dimethylsiloxane	1.0	405.1 Torr	441.1 Torr	SP	SP	SP	SP	500	1.4	1	125	20	7.7kV	0	0	0
	4-Dimethylsiloxane	1.0	405.2 Torr	441.1 Torr	SP	SP	SP	SP	400	4.25	1	100	12	7.7kV	0	0	0
	Dimethylsiloxane	1.0	405.2 Torr	441.1 Torr	SP	SP	SP	SP	125	5.1	1	100	0	7.7kV	0	0	0
	Dimethylsiloxane	1.0	405.2 Torr	441.1 Torr	SP	SP	SP	SP	5.1	1	100	100	0	7.7kV	0	0	0
	Dimethylsiloxane	1.0	405.2 Torr	441.1 Torr	SP	SP	SP	SP	5.1	1	100	80	0	7.7kV	0	0	0
	Methylsiloxane	1.0	405.2 Torr	391.1 Torr	SP	SP	SP	SP	200	5.2	1	125	10	7.7kV	0	0	0
	Methylsiloxane	1.0	405.2 Torr	391.1 Torr	SP	SP	SP	SP	5.2	1	125	20	0	7.7kV	0	0	0
	Methylsiloxane	1.0	405.2 Torr	200.1 Torr	SP	SP	SP	SP	5.2	1	100	40	0	7.7kV	0	0	0
M	Dimethylsiloxane	1.0	405.2 Torr	401.1 Torr	SP	SP	SP	SP	500	5.05	1	125	10	7.7kV	0	0	0
	Dimethylsiloxane	1.0	405.2 Torr	401.1 Torr	SP	SP	SP	SP	500	6.05	1	125	0	7.7kV	0	0	0
	Tetramethylsiloxane	1.0	405.2 Torr	811.1 Torr	SP	SP	SP	SP	6.05	1	125	10	0	7.7kV	0	0	0
	Tetramethylsiloxane	1.0	405.2 Torr	154.1 Torr	SP	SP	SP	SP	6.05	1	125	20	0	7.7kV	0	0	0
	Dimethylsiloxane	1.0	405.2 Torr	98.1 Torr	SP	SP	SP	SP	0	6.05	1	125	48	0	7.7kV	0	0
Dynamic MS/MS Parameters																	
Scan Time 3000 to										Triggered MS/MS							
										Triggered MS/MS							
										Triggered MS/MS							

Agilent's tMRM acquisition software and Veterinary Drug tMRM database and library ensure fast method implementation for targeted veterinary drug screening with compound verification.

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Match Score = 70

Ordering information:

Agilent's QQQ LC/MS Forensic Toxicology Application Kit
(G1734BA)

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